

Fever Indicator


Liquid Crystal Temperature Monitoring

The Fever Indicator is a screening tool to show preliminary results of temperature. With two indications, either a checkmark for “OK” or a sad face icon with a thermometer, you will be able to verify if someone should be concerned with a temperature above 38°C /100°F. By simply placing the Fever Indicator on a clean/dry forehead and waiting 15 seconds, color will appear next to the icon to designate temperature.

EASY-TO-READ

Through the annual flu season or a global crisis, these indicators will be able to monitor and show temperature concern over large groups, quickly and easily. With standard icons used to designate the two temperature thresholds, interpretation will be recognizable around the globe.

DIMENSIONS



C-8200 - Fever Indicator



SPECIFICATIONS

Response Temperature	Indication at 98.6°F/37°C or below Indication at 100°F/38°C or above
Size	0.625” x 0.625”
Accuracy	±1°C
Indicator Display	Reversible color change
Storage	No requirements
Adhesive	Hypoallergenic adhesive
Packaging	5,000/roll, 5 across on the roll, perf every 4 labels